The future of electrophysiology

biosignal acquisition system

General Description

The KJE-1001 multiplexed biosignal amplifier and acquisition system is a high quality, low noise level data acquisition device for 256 channels of various biological signals from freely moving animals, such as:

- electroencephalograms,
- local field potentials, and
- single neuronal activity

Due to the multiplexing capability, the information of 256 individual channels is transmitted by using only 12 wires.

The very small size of the headstage amplifiers and multiplexers allow its use on small rodents such as mice, even during freely moving behavioral tasks.



HS2 64 Channel headstage with Omnetics connectors



64 Channel headstage with onboard silicon probe



Technical Specifications

Base system:

Flexible headstage connection up to

- 256 channels
- 32 external input channels (±10 V)
- 4 digital input

Headstages:

- 0.3 Hz high pass filter
- 10 kHz low pass filter
- Gain: 400
- Input referred noise: ~1.5 μV RMS
- Sampling rate: 20 kS/Ch @ 14 bit
- Size: 14x13x3 mm, 1.9 gr (32 ch) 19x13x3 mm, 2.3 gr (64 ch) 38x34x3 mm, 6 gr (256 ch)
- +5V LED power supply

Demulitplexer:

- Time delay of demultiplexing: 50 μs
- Variable gain
- Outputs for raw signal and high pass filtered signal for MUA

System components

Base components

S-KJE-1

Amplipex KJE-1001 Basebox KJE-1001 KJE-C1A **DAQ Card** KJD-1BB **Demultiplexer Basebox Demultiplexer Unit** KJD-1000 Headstages HS3o-v1.3 32 Ch Multiplexer Headstage with Omnetics connector mates with Omnetics NPD-36-DD-GS (36 connectors + 4 guide posts) HS3h-v1.3 32 Ch Multiplexer Headstage with Hiroshe connector mates with Hiroshe DF30-FC-20DS-0.4V HS2o-v1.3 64 Ch Multiplexer Headstage with Omnetics connector mates with Omnetics NPD-36-DD-GS (36 connectors + 4 guide posts) HS2h-v1.3 64 Ch Multiplexer Headstage with Hiroshe connector mates with Hiroshe DF30-FC-20DS-0.4V HS2e-v1.3 64 Ch Multiplexer Headstage with onboard silicon probe (8x8 sites) **Cables** KJE-1cbl Shielded DAQ Interconnect Data Cable KJD-1cbl Demultiplexer Interconnect Cable HS3-cbl 32 Channel Headstage cable – Length: 12 m HS2-cbl 64 Channel Headstage cable - Length: 12 m Software



Amplipex Data Acquisition Software for 256 Channel KJE-1001 Basebox